


Protein Spotter



## “ChemPro”

A synthetic digestive agent for wet-side spotting of blood, food, eggs, milk, perspiration, vomit, discharge, and all other protein stains derived from humans or animals.

- *Effective removal of most protein stains*
- *Safe on fabrics and colors not affected by moisture*
- *Strongest protein stain remover without digestive action*
- *Flush with water or steam, level or dry before dry-cleaning*
- *Contact us for pricing and available sizes*

### DIRECTIONS

Use ChemPro<sup>®</sup> to remove blood, food, perspiration and other protein based stains. Apply ChemPro<sup>®</sup> directly to stained area. Gently work into stain with spatula or tamp with spotting brush. Flush with water or steam gun, then dry and continue the dry-cleaning process. If necessary, repeat the procedure to remove last traces of stain. Set, or extremely difficult stains may require digestive action before continuing the dry-cleaning process.

Our Family of Trusted Detergents, Spotting Agents and Specialties

